OHAUS Label Printing Instruction Manual

English

Español

Français

Deutsch

Italiano

Polski

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Operation Videos

EX Label Printing Function

OHAUS Label Designer

Table of Contents

1. INTRODUCTION	2
2. LABEL PRINTING FUNCTION	2
2.1 Step 1: Connect Explorer to label printer	2
2.2 Step 2: Enable label printing function in Explorer	3
2.3 Step 3: Print labels	3
3. PC SOFTWARE - OHAUS LABEL DESIGNER	5
3.1 Where to download the OHAUS Label Designer	5
3.2 How to use OHAUS Label Designer	5
3.2.1 Step 1: Connect Explorer to PC	5
3.2.2 Step 2: Startup	6
3.2.3 Step 3: Customize Template	7
3.2.4 Step 4: Write Label template to Explorer	11

1. INTRODUCTION

This manual is to help set up and configure an OHAUS Explorer balance (software version 2.23 and later) for printing labels and weighing results by a ZPL-compatible label printer. Please read this manual completely before software installation.

This label printing function of Explorer balance has been tested with a Zebra GC420t. It can be compatible with any label printer, which uses the ZPLII programming language.

As label formats can be highly customized, this manual is not intended as a tutorial on how to customize a label printed by the label printer.

2. LABEL PRINTING FUNCTION

Explorer series balance allows print labels in 3 steps.

- **Step 1:** Connect the balance to a label printer using an RS232 cable
- **Step 2:** Enable label printing function in Explorer
- **Step 3:** Print labels by pressing **Print** in Explorer

2.1 Step 1: Connect Explorer to label printer

For connecting an OHAUS Explorer balance to a label printer, you need to prepare the following items:


- A male-to-male RS232 cable (OHAUS PN 30057595)
- An OHAUS Explorer balance (software version 2.23 and later)
- A label printer that supports ZPL commands

Use this RS232 interface to connect the balance to a printer

Male-to-male RS232 cable

**Note:**

Usually, Explorer can be connected to a label printer with default communication settings as right.

Before printing, make sure that both balance and printer have the same settings i.e., Baud Rate, Transmission, and Handshake.

To change the connection settings of balance, enter “Menu – Communication – RS-232 Standard – Device Settings”.

2.2 Step 2: Enable label printing function in Explorer

Label Printing menu is available in Explorer balance (software version 2.23 and later) with default **Off**. To enable label printing, enter “Menu – Communication – RS-232 Standard” and set the **Label Printing** menu to **ON**.

Then **Label Template** menu appears with default English label template “100mm * 75mm”.

2.3 Step 3: Print labels

To print labels, choose label template and press **Print** button in Explorer.

There are 2 built-in templates in Explorer.

- **#1 - 100mm * 75mm**
- **#2 - 100mm * 10mm**

Below are **printout examples** using default templates:

Template	#1 - 100mm * 75mm	#2 - 100mm * 10mm
Size	<p>100 mm</p> <p>75 mm</p>	<p>100 mm</p> <p>10 mm</p>
Content	<ul style="list-style-type: none"> • Data/Time • Sample ID • Sample Name • Batch ID • Result • Gross Weight • Tare Weight • Net Weight • User Name 	<ul style="list-style-type: none"> • Result

When Explorer is in label printing mode, it prints stable results only.

Please refer to **Section 3 PC Software - OHAUS Label Designer** for instruction to change the default label templates.

3. PC SOFTWARE - OHAUS LABEL DESIGNER

OHAUS Label Designer enables the user to edit label templates and save to Explorer by connecting to a computer.

3.1 Where to download the OHAUS Label Designer

Click to get free [OHAUS Label Designer](#) software.

You can also download the latest software from the official website **OHAUS.com** at the path:

3.2 How to use OHAUS Label Designer

3.2.1 Step 1: Connect Explorer to PC

This software works with the software version of Explorer 2.23 and later.

The operating system of PC for running OHAUS Label Designer is Windows 10® 32 / 64 bit.

To connect the OHAUS Explorer balance with PC, you need the below accessory.

Interface kit RS232-USB (OHAUS PN 30304101)

3.2.2 Step 2: Startup


Figure 3.2.2-1

1. Ensure the balance has been logged in as Admin or Supervisor before running the OHAUS Label Designer software.
2. Double-click the software.
3. Choose the operation Language in the upper right corner. Below languages are available.
 - English
 - Spanish
 - French
 - German
 - Italian
 - Russian
 - Japanese
 - Chinese
4. Check the communication settings.
 - Port number is automatically detected after connection.
 - The PC software and balance have the same default communication Setting (see figure 3.2.2-1).
5. Click **OK** to continue.
6. Once the balance has successfully connected to the PC, it shows the balance **Model** and **Serial Number** (see figure 3.2.2-2).



Figure 3.2.2-2

3.2.3 Step 3: Customize Template

Main screen

After startup, OHAUS Label Designer appears with a blank state in the editing zone. The main screen has 4 sections.

1
Import label template from PC /
Export label template to PC
in *.tmp format

2
Available values (print contents)
read from Explorer


3
Zone for editing label template
or showing imported template

4
Read label template from
Explorer / Write label template
to Explorer


Customize label format

All label templates should be edited in ZPL programming language. **Maximum 100 lines with up to 80 characters in each line can be used in each label template.**

To create a new label template, follow these steps:

1. Click **Read**, select Print Template, for example, #1 – 100mm * 75mm. Then click **OK** to continue.

2. The below items can be customized by ZPL scripts
 - text position
 - font size
 - customized contents
 - Values read from balance (in blue text)

Example 1:

To add “**date**” and “**time**” into the label, click **[5200] Date** and **[5201] Time** from the left side; and on the right side, ZPL scripts will show the code in blue text `[$5200!] [$5201!]` at the line.

Printout	ZPL Scripts
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> Date/Time: 12/12/2021 05:39:23 </div>	<pre>^FT130,68^A0N,28,28^FH^FDDate/Time:^FS ^FT447,68^A0N,28,28^FH^FD[\$5200!] [\$5201!]^FS</pre>

You can also define the font size of “date and time” and its position in the label.

- **28, 28** indicate the font size (height and width) of both “Date/Time” and `[$5200!] [$5201!]` (values read from Explorer).
- **Date/Time:** can be customized by users. For example, you can change it to **Date and Time:**.
- **130,68** indicate the position of **Date/Time:** in the label, i.e, from the left to right 130 dots, from top to down 68 dots. According to the label size, you can adjust the position, i.e, **150, 88**.

ZPL Scripts	<code>^FT130,68^A0N,28,28^FH^FDDate/Time:^FS</code>	<code>^FT150,88^A0N,28,28^FH^FDDate and Time:^FS</code>
Printout	<p>Diagram illustrating the printout of "Date/Time:". The text is positioned at 130 dots from the left and 68 dots from the top.</p>	<p>Diagram illustrating the printout of "Date and Time:". The text is positioned at 150 dots from the left and 88 dots from the top.</p>



Example 2:

To print "**Sample ID**" as barcode, use below ZPL scripts and click **[5311] Sample ID** from left side, and on the right side, ZPL scripts will show the code at the line of "`^BY2,3,36^FT,447,129^BCN,,Y,N ^FD>${5311}!^FS`".

Printout	ZPL Scripts
	<code>^BY2,3,36^FT,447,129^BCN,,Y,N ^FD>\${5311}!^FS</code>

You can define the text position of bar code in the label.

- `447,129` indicate the poision of `${5311}!` (barcode) in the label, i.e, from the left to right 447 dots, from top to down 129 dots. According to the label size, you can adjust the position, i.e, `130, 129`.

ZPL Scripts	<code>^BY2,3,36^FT,447,129^BCN,,Y,N ^FD>\${5311}!^FS</code>	<code>^BY2,3,36^FT,130,129^BCN,,Y,N ^FD>\${5311}!^FS</code>
Printout		

To print "**Sample ID**" as number, use below ZPL scripts and click **[5311] Sample ID** from the left side, and on the right side, ZPL scripts will show the code at the line of "`^FT447,129^A0N,28,28^FH^FD${5311}!^FS`".

You can define the font size of sample ID and the text position (see **Example 1**).

Printout	ZPL Scripts
	<code>^FT447,129^A0N,28,28^FH^FD\${5311}!^FS</code>

Reset

You can reset the edited label template to default template (#1 - 100mm*75mm) by clicking **Reset**.

Model:


Serial Number:

Print Template:

3.2.4 Step 4: Write Label template to Explorer


5 label templates (#1 - #5) are allowed to store in Explorer balances.

To write templates to the balance, follow below steps:

- Click **Write**.
- Choose the template # you want to save to and define the template name. Click **OK** to continue.


For example, if you choose #3, the original format in #3 will be replaced by new format automatically.

- The left screen will appear if the template is written to Explorer successfully.

To backup more templates to your PC, click **Export** and save the templates in format (*.tmp).

These templates can be imported to OHAUS Label Designer for editing at any time.

For additional support for Zebra printers and label formats, contact Zebra directly at www.zebra.com.

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P/N 30708369 B © 2021 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés /
Ohaus Corporation